



MEDICINE REMINDER AND PILL DISPENSER

GROUP 07

Our Team

CHATHUMINI V.G.A.

240078J

GURULLAWALA G.G.Y.M.

240210J

KASTHURIARACHCHI I.D.

240338H

NAVEENDRAN S.

240453F

SARAWANADASAN P.P.

240599G

Introduction

- Patients with dementia often experience memory loss and difficulty following medication schedules.
- Missing or incorrect medication doses can lead to serious health complications.
- Caregivers need an efficient way to monitor the regular medicine use of patient.

What Makes Our System Different?

- Automatically dispenses the exact prescribed number of pills using Archimedes' screw mechanism.
- Supports up to five different medicine types within a single system.
- Verifies pill dispensing using a VL53L0X sensor for improved dosage accuracy.
- Provides automatic water dispensing together with medication.
- Detects medicine collection using an IR sensor
- Enables real-time caregiver monitoring through a mobile application.

Our system integrates medication reminders, pill dispensing, water supply, and caregiver monitoring into a single smart healthcare solution.




Project Goal and Objectives

Goal: Develop a smart medicine reminder and automated pill dispensing system that improves medication adherence among dementia patients.

- Objectives:**
- Remind patients at scheduled medication times.
 - Dispense the correct number of pills automatically.
 - Supply drinking water together with the medication.
 - Monitor medicine collection by the patients.
 - Notify caregivers through a mobile application.
 - Reduce medication errors and missed doses.

Our System

- Stores five different medicine types.
- Dispenses the required number of pills automatically.
- Provides medicine reminders through an alarm.
- Supplies water automatically.
- Monitors medicine collection.
- Sends notifications to caregivers.



Hardware Components



ESP32




12V POWER ADAPTER



IR SENSOR



ULTRASONIC SENSOR



VL53L0X SENSOR



NEMA17
STEPPER MOTOR



28BYJ-48
STEPPER MOTOR



SERVO MOTOR




PUMP

Software Tools



ARDUINO IDE




ANDROID STUDIO



SOLIDWORKS

How It Works?






Pill Dispensing Mechanism

Pill Dispensing

- Each medicine stored in a separate container.
- The container uses an Archimedes' screw driven by a stepper motor.
- The number of motor rotations determines the dosage.

Conveyor Belt

- The pills move into the medicine cup after dispensing.
- VL53L0X counts the number of pills passing through it.
- A servo motor removes extra pills to continue the medicine supply.




Monitoring and Notification System

➤ App-controlled dosage and notifications.

Notifications are sent to the caregiver:

- When medication time reaches the scheduled value.
- When it confirms the medicine collection by the IR sensor.
- When tracking the medicine content in each container after each dosage.
- When the water level of the tank reaches the minimum value.



Water Dispensing System

Components

- Water tank
- Water pump
- Ultrasonic sensor

Functions

- Dispense sufficient water automatically.
- Monitor the water level.
- Alert caregivers when the water level is low.

Advantages

- Improves medication adherence.
- Reduces missed medication doses.
- Supports dementia patients.
- Provides accurate dosage dispensing.
- Allows caregiver monitoring.
- Reduces human errors.
- Increases patients' independence.


Limitations

- A limited number of pills can be contained in each spiral.
- Supports only medicines compatible with the dispensing mechanism.
- Requires a continuous power supply.
- Limited to five medicine containers.
- Cannot directly verify medicine consumption after collection.

Future Improvements

- A pill sorting system to automatically load pills into the spirals.
- An efficient safety mechanism using image recognition.
- Expanded pill containers for more doses.
- Integration with hospital databases.
- Voice-assisted reminders.
- Battery backup system.

Conclusion



Successfully
Developed



Reminders And
Monitoring



Improved Patient
Care



Thank You!